

SCREWING ACTIONS FOR EACH PLACED IMPLANT

		Omplant DRIVER or hand driver (18 Ncm)										RATCHET TORQUE (30 Ncm)													
	*Surgical Phase			**Prosthetic Phase																					
				Impression			Metal fit Check				Porcelain fit Check				Delivery										
Implant	1 IN cover screw	2 OUT cover screw	3 IN healing abutment	4 OUT h.a	5 IN organizer	6 IN transfer screw	7 OUT	8 OUT organizer	9 IN h.a	10 OUT h.a	11 IN	13 IN	14 OUT	16 OUT organizer	17 IN h.a	18 OUT h.a	19 IN	21 IN	22 OUT	24 OUT organizer	25 IN h.a	26 OUT h.a	27 IN	29 IN	30 IN
Analog							abutment screw				abutment screw				abutment screw										

*More clinical usages:

1. screwing cover screws or healing abutments to unseen implants (partially covered with gums)
2. locating failed implant - (not osseointegrated)

** NOT included : false in/out trials!, metal welding need, or changing healing abutment during the transitional restoration